

WHEELCHAIR LEGRESTS/FOOTRESTS

HEAVY DUTY ARTICULATING LEGRESTS

Everest & Jennings introduces the heavy duty articulating elevating legrests, industry compatible with aluminum footplate for durability. Normally used by tall individuals who require extension while fully elevated, this elevating legrest is sold as a separate accessory.

Features

- Aluminum footplate in black for bariatric or standard manual wheelchairs
- Extends in length from a minimum 16" in the downward position up to a maximum of 27" when fully elevated
- Swingaway detachable legrest can be adjusted for different needs or to avoid fatigue and enhance treatment
- Chrome-plated finish with black padded panel in vinyl for comfortable calf support
- Compatible with Advantage, Advantage Recliner, Traveler SE, Traveler LX, Traveler L3, Traveler L3 Plus, Traveler L4, Traveler HD & Traveler HTC
- Meets Medicare code K0053

Warranty

One (1) year limited warranty



Length Adjustment Range*

	Position Downward min-max		Fully Elevated Upward min-max		UOM
	inch	mm	inch	mm	
90763432	16-22	40.6-55.9	22-27	55.9-68.6	1 pr

*Measurements from seat upholstery edge to footplate



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2016, GF Health Products, Inc. All Rights Reserved.
Everest & Jennings is a registered trademark of GF Health Products, Inc.

graham
field

GF1600124RevB55