

FR601E

LUMEX®
by graham-field

The Powered Pivot-Arm Clinical Care Recliner our flagship recliner, is designed for patient comfort and caregiver convenience. The FR601E offers powered infinite recline positions, pivoting arms, and industry-leading comfort. With high-performance viscoelastic memory foam in the seat and backrest, it delivers comfort, outstanding support, and pressure redistribution. Recliner has Trendelenburg positioning, a manual recline override, and battery backup for operation during power failure. Optional features include USB charging ports (Type A & C), Central Locking Casters, and Foley Hook.



CLINICAL
CARE



**500
POUND™**
WEIGHT CAPACITY
EVENLY DISTRIBUTED

PACU	Rehab	Oncology	IV Therapy	Exam Room	Patient Room	Emergency Room	Dialysis (Facility)	Dialysis (In Home)	Nocturnal Dialysis	Physical Therapy	Pain Management	Post Operative Care	Same Day Surgery / ASC
☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒

AREAS OF USE

HEALTHCARE SEATING SOLUTIONS

**graham
field™**

Options:



FULLY UPHOLSTERED ARMS
FR601EPxxxx-xxxx



URETHANE ARM CAPS
FR601EUxxxx-xxxx



HEAT ONLY
FR601Exxx0xx-xxxx



HEAT AND MASSAGE
FR601ExxxHxx-xxxx



FOLEY HOOK
FR601ExxxxxF-xxxx



CENTRAL LOCKING CASTERS
FR601ExCxxx-xxxx



USB-A + USB-C CHARGE PORTS
FR601ExxxBx-xxxx

† -xxxx in the item number is the placeholder for upholstery code

FR601E Product Info



Positions: POWERED INFINITE POSITION + CAREGIVER-ACTIVATED TRENDLENBURG

Infinite Power Recline Positions: Allows Patient to Set a Position That is Most Comfortable
Trendelenburg Position: Caregiver-Activated Using Rear Mounted Hand Trigger



Features:

- Both arms pivot outward, enabling easier cleaning and access
- Aesthetically pleasing design
- Wide contoured seat and back provides an enhanced level of patient comfort
- Viscoelastic Memory Foam engineered into the seat and backrest for ultimate comfort
- 5-inch Tente® swivel casters for ease of maneuverability (also available with Central Locking Casters)
- Two self-storing, fold down, side tables with smooth surface for easy cleaning
- Comfortable Urethane armrest caps are easy to maintain (also available with fully upholstered arms)
- Pullout footplate aids in patient entry and exit and safer patient transport
- Rear push handle for improved patient transport



FOLD DOWN, SIDE TABLES



PULLOUT FOOTPLATE



TENTE® CASTERS



PIVOT-ARM RELEASE



REAR PUSH HANDLE



TRENDELBURG TRIGGER



HAND REMOTE

Accessories:

★ sold separately	item no. [†]
Headrest Pillow	FR5659xxxUS
Headrest Skirt	HRSxxx
Headrest Cover	HRC601xxx
IV Pole	2101A
IV Mount ⁵	FR601IVMNT

Specifications:

Overall Width						Overall Depth		Overall Height		Arm Width		Width between arms		Arm Height from floor		Seat Width		Seat Depth		Seat Height		Back Width		Back Height above seat		Weight ¹ Capacity			
tables down		central locking		side tables up																									
inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	lb	kg†		
33.00	83.82	35.00	88.90	53.50	135.89	42.00	106.68	72.00	182.88	51.00	129.54	3.00	7.62	24.00	60.96	28.50	72.39	23.00	58.42	19.00	48.26	24.00	60.96	30.50	77.47	30.00	76.20	500	226.80

*Ships Freight: request additional shipping quote

† -xxxx in the item number is the placeholder for upholstery code

† standard rounding method used

¹ maximum weight capacity - evenly distributed - refer to user manual

⁵ requires factory installation include item with purchase

★ options available at additional cost - accessories sold separately

Available Colors:

COLOR OPTIONS AND FABRIC SWATCHES ON WEBPAGE Upholstery Selections



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Images are digital reproductions. perceived texture and coloring may vary on finished product. Monitor and printer settings will effect appearance of color.

• 1 770.368.4700

Information contained herein is subject to change without notice. The most current & complete product information can be found on our website.

www.grahamfield.com



© 2026, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.
Tente is a registered trademark of TENTE International GmbH.

graham field

GF2600107/Rev/A26