



 **SureStryde™**

*Elevate Your Journey and
Enjoy Life's Richest Moments*



LUMEX®
by graham-field

Features

- Patented seat slides along frame allowing user to walk inside the rollator frame allowing for improved gait
- Adjustable push bar handle provides a place to propel rollator forward with better posture and stability
- Ergonomic handgrips with easy brake access
- Intuitive one brake handle can activate both brakes
- Contoured, padded backrest offers comfort when seated
- Deluxe, soft pouch for personal items can be easily installed and removed
- Sleek frame design makes rollator functional and versatile to meet various user needs
- Large wheels allow for easy maneuverability
- Maximum weight capacity 300 lb (136 kg) **EVENLY DISTRIBUTED**
- Folding mechanism allows rollator to fold for storage and transport
- Bonus curb climber included



	overall Width		overall Depth		Width inside handgrips		minimum Handle Height		maximum Handle Height		Wheel front		Wheel rear	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
S9000S	24.5	62.2	31.7	80.5	16	40.6	29.4	74.7	35.4	89.9	8	20.3	9.5	24.1

	Seat Width		Seat Depth		Seat Height from floor		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM	Color	Curb Climber
	inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg			
S9000S	17.2	43.7	15.1	38.4	22	55.9	22	9.98	300	136	1/ea	Silver Storm	Yes



Bright Reflectors



Comfortable Handgrips



Curb Climber



Collapsible



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

***1 770.368.4700**

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website
www.grahamfield.com



© 2024, GF Health Products, Inc. All Rights Reserved.
Lumex is a trademark of GF Health Products, Inc.

grahamfield

GF2400117RevA24